Kit de HBcAb en ELISA

Cat. No. ITP21501

USO PREVISTO

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BIBLIOGRAFIA

1. Engvall E, and Perimann P. Enzyme-linked immunosorbent assay (ELISA) Quantitative assay of immunoglobulin G. *Immunochem*, 8:871-874, 1971
2. Engvall E, and Perimann P. Enzyme-linked immunosorbent assay (ELISA) Protides of the Biological Fluids, Proceedings of the Nineteenth Colloquium. Brugge (Peeters H., ed.) Pergamon Pre, Oxford, PP. 553-556, 1971.
3. Engvall E, Jonsson K, and Perimann P. Enzyme-linked immunosorbent Assay II. Quantitation Assay of Protein antigen, immunoglobulin G. by means of enzyme-labeled antigen and antibody-coated tubes. *Biochim. Biophys. Acta* 251: 427-434, 1971.
4. VanWeemen BK, and Schuurs AHWM Immunoassay using antigen-enzyme conjugates, *FEBS Letters* 15: 232-236, 1971.
5. Wisdom GB. Enzyme-immunoassay, *Clin. Chem.* 22: 1234-1255, 1976.
6. Wolters G, Kuijpers, L, Kacaki J, and Schuurs A. Solid-phase enzyme-immunoassay for detection of hepatitis B surface antigen, *J. Clin. Pathol.* 29:873-879. 1976.
7. Wei R, Knight GJ, Zimmeman DH, and Bond HE, Solid-phase enzyme-immunoassay for detection of hepatitis B surface antigen, *Clin. Chem.*, 23:813-815, 1977.
8. David GS, Present W, Martinis J, Wang R, Bartholomew R, Desmond W, and Sevier ED, Monoclonal antibodies in the detection of hepatitis infection, *Med. Lab. Sci.* 38:341-348, 1981.
9. Drouet J, Courouze A-M, Kalil J., Ropars C, Felious M, Monoclonal antibodies to HBsAg produced by murine hybridomas. (in) *Viral Hepatitis* (Szmuness, W., Alter, H.J., Maynard, J. E., eds.) Franklin Institute Press, Philadelphia, pp. 706-707, 1982.
10. Goodall AH, Miescher G, Meek FM, Janossy G, Thomas G, Thomas HC. Monoclonal antibodies in a solid-phase radiometric assay for HBsAg *Med. Lab. Sci.* 38: 349-354, 1981.
11. Kennedy RC, Loniscu-Matiu I, Alder-Alder-Storthz K, Henkel RD, Sanchez Y, Dreesman GR. Characterization of Anti-Hepatitis B Surface Antigen Monoclonal Antibodies, *Intervirology.* 19: 176-180, 1983.
12. Shih JW-K, Cote PJ, Dapolito GM, and Gerin JL. Production of monoclonal antibodies against Hepatitis B surface antigen (HBsAg) by somatic cell hybrids, *J. Virol. Meth.* 1: 257-273, 1980.

13. Wands JR, Zurawski VR, High Affinity Monoclonal Antibodies to Hepatitis B Surface Antigen (HBsAg) Produced by Somatic Cell Hybrids, *Gastroenterology* 80: 225-232, 1981.
14. US Environmental Protection Agency. EPA guide for infectious waste management. Washington, DC: May, 1986.
15. Centers for Disease Control Recommendations for prevention of HIV transmission in health care settings. *MMWR* 36, Supplement No. 2S, 1987.
16. Sehulster LM, Hollinger FB, Dreesman GR, and Melnick JL. Immunological and biophysical alteration of Hepatitis B virus antigens by sodium hypochlorite disinfection. *Appl. and Envir. Microbial.*, 42: 762-767, 1981.
17. Mehan BH, University Chemistry, Second Edition, Addison Wesley, p. 648, 1968.
18. Grangeot-Keros L, Lamber T. Dubreuil P, Briantais MJ, and Pillot J: False Reactions in Radioimmunoassay for Viral Hepatitis B Markers in Patients Suffering from Coagulation Disorders. *Vox. Sang:* 42: 160-163, 1982.
19. Technical Manual of the American Association of Blood Banks, Eighth Edition, Philadelphia, J.B., Lippincott Co., p. 216, 1981.
20. Epidemiologic Notes and Reports, Hepatitis B Contamination in a Clinical Laboratory-Colorado, *MMWR* 29: 459-565, 1980.

InTec PRODUCTS, INC.