



ENDLESS CHALLENGE
ENDLESS IMPROVEMENT

BA 680

Auto Urine Analyzer



BA 680

Auto Urine Analyzer



- **Test items** : 11 parameters : Glucose, Bilirubin, Ketone, Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite, Leukocyte, Vitamin C
- **Test method** : Dry strip color-developing method
- **Test principle** : Dual wavelength reflectance photometry
- **Light source** : Cold lamp with high brightness and long lifetime
- **Test system** : Color-indentify system with high sensitivity and repeatability
- **Wavelength** : 520nm, 640nm, 660nm
- **Test mode** : Single and continuous test mode are selectable
- **Test slot** : Using sliding slot, easy to clean and replace, low carryover
- **Throughput** : 500 Strips/Hour
- **Waste collection** : Waste urine and waste strips can be collected automatically
- **Test strip** : Use GENIUS original strips only.
- **Display** : LCD with high brightness and resolution
- **Language** : Customizable language software
- **Memory** : Save 50000 samples results, and auto recover the results when power off during working
- **Interface** : RS-232 serial port, can be connected to external printer and computer
- **Printer** : Inner thermal printer, can be connected to external printer
- **Report format** : Half quantitative and qualitative, report the color and turbidity is available
- **Result unit** : SI unit and routine unit are selectable
- **Status check** : Auto-check the status during power on and off, auto-indentify the alarm
- **Working environment** : Temperature: 10°C ~ 30°C; Humidity: 20%RH ~ 80%RH
- **Voltage** : 85-264V, Auto adapt according to the input voltage
- **Power** : ≤80W
- **Quality system** : Recognized by TUV and CMD

